

273 | Test Tube Aging Tester Model TA-1



■ APPLICATION

Metal bath for conducting heat deterioration or aging test by keeping the sample inside a glass tube for eliminating the effect of migration of additives, etc.

■ SPECIFICATIONS

Model	TA-1
Max. temperature	300°C
Test tube hangings	24 (test tubes are supplied as an option) (Diameter 38mm x Length 300mm) <i>Note: Glass tubes are provided as optional</i>
Test tube dimensions	Ø38 x h300mm
Power requirement	Single-phase, AC200V, 50Hz or 60Hz, 5.0kVA
Dimensions	W500 x D720 x H450mm
Weight (approx.)	120kg
Related standards	ISO 188 JIS K 6257 ASTM D865

■ OPTIONS

Model	Name	Remarks
TA-GT	Test tube set	Set of following accessories: Cork (24pcs.), specimen hook (24pcs.), test tube (24 pcs.), washer (50pcs), aluminum tape
PA-1	Heater ON/OFF & buzzer OFF option	

Specifications are subject to change without notice.



TOYO SEIKI SEISAKU-SHO, LTD.

5-15-4, Takinogawa, Kita-ku, Tokyo 114-8557, Japan

Tel: +81-3-3916-8183 Fax: +81-3-3916-8173

www.toyoseiki.co.jp

20250317 MN (2024.09 JP)